

RESISTANCE PROFILE, TERBINAFAINE
RESISTANCE SCREENING AND MALDI-TOF MS
IDENTIFICATION OF THE EMERGING PATHOGEN

TRICHOPHYTON INDOTINEAE

Ir Roelke De Paepe

Mycology & Aerobiology unit, Sciensano, Brussels, Belgium

Introduction

- Epidemic of tinea infections in (North) India – mainly *tinea corporis* and *tinea cruris*
- Patients demonstrate highly inflammatory, severe lesions
- High level of resistance to antifungal agent **terbinafine** (TRB)
- Emerging prevalence in neighbouring countries, as well as world-wide



Uhrlaß et al., 2022. *Trichophyton indotineae*-An Emerging Pathogen Causing Recalcitrant Dermatophytoses in India and Worldwide-A Multidimensional Perspective

***Trichophyton mentagrophytes* ITS genotype VIII** was identified as causal agent of these infections

Introduction

T. mentagrophytes ITS genotype VIII or *T. indotineae*?

- Morphologically identical to *T. mentagrophytes*, genetically almost identical to *T. interdigitale*
- But: phylogenetic clustering of the 'Indian' isolates
ecological, physiological and clinical differences
high resistance level to TRB

Kano *et al.*, 2020:
new species ***T. indotineae***

→ **Correct diagnosis is often difficult**, especially when only using morphology for identification


Introduction

High level of **terbinafine resistance** within *T. indotineae* strains

- Consequence of excessive use of over-the-counter antifungal creams and topical steroids or glucocorticoids?
- Point mutations in the **squalene epoxidase-encoding gene (SQLE)**
 - squalene epoxidase (SQLE) is an essential enzyme in the biosynthesis of ergosterol
 - most common amino acid substitutions: **F397L** and **L393F**
- Resistance to other antifungal agents as well?



Research




Terbinafine resistance screening




①
Testing a **TBF containing screening medium** as a quick way to determine TBF resistance of *T. indotineae* strains

4 DAYS:




7 DAYS:



- Screening on Sabouraud medium amended with TRB (0.2 µg/mL)
- Incubation of *T. indotineae* strains (n=20) at 25°C
- After **4 days** of incubation: fungal growth of **resistant strains**, while no visible growth for sensitive strains

Figure 1. Growth of the *T. indotineae* isolates (n=20) on a screening medium for terbinafine resistance. The medium consists of Sabouraud dextrose agar, amended with terbinafine.

Antifungal susceptibility testing



②

Performing **antifungal susceptibility testing** (AST, following the EUCAST method) to check for possible multi-resistance

- EUCAST broth microdilution antifungal susceptibility testing
- Minimal inhibitory concentrations at 100% inhibition of fungal growth (MIC100):

fluconazole	terbinafine
itraconazole	griseofulvin
voriconazole	ciclopirox olamine
ketoconazole	naftifine
	amorolfine




EUCAST

EUROPEAN COMMITTEE
ON ANTIMICROBIAL
SUSCEPTIBILITY TESTING

European Society of Clinical Microbiology and Infectious Diseases

- Tested concentrations ranged from 0.008 – 64 µg/mL
- WT-UL values for itraconazole, voriconazole, terbinafine and amorolfine
TRB: 0.25 µg/mL, ITC: 1.0 µg/mL, VOR: 2.0 µg/mL and AMO: 1.0 µg/mL

Antifungal susceptibility testing



②

Performing **antifungal susceptibility testing** (AST, following the EUCAST method) to check for possible multi-resistance

Table 1. MIC-values ($\mu\text{g/mL}$) of all *T. indotineae* strains ($n = 20$) for tested antifungal agents fluconazole (FLC), itraconazole (ITC), voriconazole (VOR), ketoconazole (KET), terbinafine (TRB), griseofulvin (GRI), ciclopirox olamine (CPX), naftifine (NAF) and amorolfine (AMO).


<i>T. indotineae</i> strains	TRB resistance	MIC ($\mu\text{g/mL}$)								
		FLC	ITC	VOR	KET	TRB	GRI	CPX	NAF	AMO
IHEM 28378	yes	16	0.032	0.25	0.25	16	8	1	32	0.25
IHEM 28379	yes	32	0.125	0.25	0.25	16	2	1	32	0.125
IHEM 28380	yes	64	0.25	0.25	0.5	16	4	1	32	0.125
IHEM 28384	yes	16	0.125	0.25	0.25	16	4	1	32	0.125
IHEM 28386	yes	32	0.125	0.25	0.25	16	4	0.5	32	0.125
IHEM 28387	yes	64	0.25	0.5	1	16	8	0.5	32	0.25
IHEM 28392	yes	64	0.25	0.25	0.25	16	4	0.5	32	0.125
IHEM 28394	yes	32	0.125	0.25	0.125	16	2	0.5	16	0.25
IHEM 28395	yes	32	0.125	0.25	0.25	16	4	0.5	32	0.125
IHEM 28381	no	64	0.25	0.25	0.5	0.016	2	1	0.25	0.032
IHEM 28382	no	32	0.125	0.125	0.5	0.064	2	1	0.125	0.064
IHEM 28383	no	64	0.25	1	1	0.064	4	1	0.25	0.064
IHEM 28385	no	32	0.5	0.5	0.5	0.032	4	0.5	0.125	0.064
IHEM 28388	no	8	0.125	0.125	0.125	0.064	2	1	0.064	0.064
IHEM 28389	no	16	0.125	0.125	0.125	0.125	2	1	0.125	0.125
IHEM 28390	no	64	0.25	0.5	0.5	0.125	4	1	0.25	0.125
IHEM 28391	no	64	0.25	0.5	0.5	0.032	4	0.5	0.25	0.064
IHEM 28393	no	16	0.064	0.125	0.125	0.064	8	0.5	0.5	0.25
IHEM 28396	no	32	0.25	0.25	0.5	0.032	2	0.5	0.25	0.25
IHEM 28397	no	16	0.064	0.125	0.125	0.032	2	0.5	0.064	0.125



Table 2. Geometric means (GM) of the MIC-values of all *T. indotineae* strains, and separately for the terbinafine (TRB) resistant and susceptible strains. Asterisks (*) symbolize a significant difference between the GM MIC-values of TRB susceptible and resistant strains, after statistical analysis (Mann-Whitney-Wilcoxon test, $p = 0.05$).

Antifungal agent	GM MIC-values ($\mu\text{g/mL}$)		
	All ($n = 20$)	TRB susceptible ($n = 11$)	TRB resistant ($n = 9$)
Fluconazole	33.13	30.05	37.33
Itraconazole	0.16	0.17	0.15
Voriconazole	0.25	0.25	0.25
Ketoconazole	0.31	0.32	0.29
Terbinafine	0.67	0.05*	> 16.00*
Griseofulvin	3.48	2.92	4.32
Ciclopirox olamine	0.71	0.73	0.68
Naftifine	1.75	0.17*	29.63*
Amorolfine	0.12	0.09*	0.17*

Antifungal susceptibility testing



②


Performing **antifungal susceptibility testing** (AST, following the EUCAST method) to check for possible multi-resistance

Antifungal agent	GM MIC-values (µg/mL)		
	All (n = 20)	TRB susceptible (n = 11)	TRB resistant (n = 9)
Fluconazole	33.13	30.05	37.33
Itraconazole	0.16	0.17	0.15
Voriconazole	0.25	0.25	0.25
Ketoconazole	0.31	0.32	0.29
Terbinafine	0.67	0.05*	> 16.00*
Griseofulvin	3.48	2.92	4.32
Ciclopirox olamine	0.71	0.73	0.68
Naftifine	1.75	0.17*	29.63*
Amorolfine	0.12	0.09*	0.17*

Table 2. Geometric means (GM) of the MIC-values of all *T. indotineae* strains, and separately for the terbinafine (TRB) resistant and susceptible strains. Asterisks (*) symbolize a significant difference between the GM MIC-values of TRB susceptible and resistant strains, after statistical analysis (Mann-Whitney-Wilcoxon test, $p = 0.05$).

- Significant increase in geometric mean of the minimal inhibitory concentration (MIC) between TRB susceptible and resistant strains for **naftifine** and **amorolfine**
- Naftifine and terbinafine are both allylamines, while amorolfine is a morpholine antifungal agent
- All strains were considered susceptible for itraconazole, voriconazole and amorolfine using WT-UL values by EUCAST
- No significantly elevated MIC values for azoles


DNA sequencing of *SQLE*



MIC (µg/mL)												SQLE missense mutations	amino acid substitution
<i>T. indotineae</i>	TRB		ITC	VOR									
known to cause TRB resistance			0.032	0	known to cause azole resistance							T1189C	F397L
IHEM 28379	yes	32	0.125	0.25	0.25	16	2	1	32	0.125		T1189C	F397L
F397L L393F TBF resistance & possible cross-resistance for naftifine A448T Not specifically leading to TBF or azole resistance on its own												T1189C	F397L
												T1189C	F397L
												T1189C	F397L
												C1191A	F397L
												T1189C	F397L
												C1191A	F397L
												A1179C	L393F
												T1189C & G1342A	F397L & A448T
												G1342A	A448T
												none	none
												G1342A	A448T
												G1342A	A448T
												none	none
												G1342A	A448T
												G1342A	A448T
												G1342A	A448T
IHEM 28396	no	32	0.25	0.25	0.5	0.032	2	0.5	0.25	0.25		G1342A	A448T
IHEM 28397	no	16	0.064	0.125	0.125	0.032	2	0.5	0.064	0.125		G1342A & G1330A	A448T

Table 1. MIC-values (µg/mL) of all *T. indotineae* strains (n = 20) for tested antifungal agents fluconazole (FLC), itraconazole (ITC), voriconazole (VOR), ketoconazole (KET), terbinafine (TRB), griseofulvin (GRI), ciclopirox olamine (CPX), naftifine (NAF) and amorolfine (AMO). MIC-values were obtained by following EUCAST guidelines for the antifungal susceptibility testing of dermatophytes. Missense mutations in *SQLE* and the corresponding amino acid substitutions in the *SQLE* protein for the tested TRB resistant and susceptible strains. Nucleobases: thymine (T), cytosine (A), adenine (A), guanine (G). Amino acids: phenylalanine (F), leucine (L), alanine (A), threonine (T).

Evaluation of MALDI-TOF MS as an identification tool



- Matrix-Assisted Laser Desorption/Ionization Time Of Flight Mass Spectrometry
- Protein extraction based on protocol described by Cassagne *et al.* (2011)
- Spectra were compared to two databases:

① in-house database BCCM/IHEM Sciensano



② MSI V2.0 database
<https://msi.happy-dev.fr/>



- Can MALDI-TOF MS differentiate TRB resistant and susceptible strains?
Do susceptible strains have a susceptible reference strain as a best-match and vice versa?

Evaluation of MALDI-TOF MS as an identification tool





Figure 2. MALDI-TOF MS identification percentages on species and genus level for both databases.

- MALDI-TOF MS can **reliably distinguish** *T. indotineae* strains from closely related species, using both an in-house and public spectral database.
- It could **not differentiate** between TRB susceptible and resistant *T. indotineae* strains.




Figure 3. Column charts of the best-matched reference strains of the *T. indotineae* strains.

Conclusions

- A terbinafine containing Sabouraud medium could screen *T. indotineae* strains for TRB resistance in just four days.
- TRB resistant *T. indotineae* strains showed significantly elevated MIC values for naftifine and amorolfine, raising suspicions of cross-resistance (for naftifine).
- The observed TRB resistance was caused by point mutations in *SOLE*, resulting in amino acid substitutions F397L and L393F, while amino acid substitution A448T does not seem to cause TRB nor azole resistance when no other substitutions are present.
- MALDI-TOF MS is a reliable tool for the identification of *T. indotineae* strains, as it is able to distinguish them from closely related species *T. mentagrophytes* and *T. interdigitale*.
- MALDI-TOF MS could not distinguish TRB resistant *T. indotineae* strains from susceptible ones.

Acknowledgments



Dr Pietro Nenoff
Dr Silke Urhlaß



Anne-Cecile Normand
Dr Renaud Piarroux

References

- Uhrlaß, S., Verma, S. B., Gräser, Y., Rezaei-Matehkolaei, A., Hatami, M., Schaller, M., & Nenoff, P. (2022). *Trichophyton indotineae*-An Emerging Pathogen Causing Recalcitrant Dermatophytoses in India and Worldwide-A Multidimensional Perspective. *Journal of fungi (Basel, Switzerland)*, 8(7), 757. <https://doi.org/10.3390/jof8070757>
- Kano, R., Kimura, U., Kakurai, M. *et al.* *Trichophyton indotineae* sp. nov.: A New Highly Terbinafine-Resistant Anthropophilic Dermatophyte Species. *Mycopathologia* 185, 947–958 (2020). <https://doi.org/10.1007/s11046-020-00455-8>
- Cassagne C, Ranque S, Normand AC, Fourquet P, Thiebault S, et al. (2011) Mould Routine Identification in the Clinical Laboratory by Matrix-Assisted Laser Desorption Ionization Time-Of-Flight Mass Spectrometry. *PLOS ONE* 6(12): e28425. <https://doi.org/10.1371/journal.pone.0028425>
- Yamada T, Maeda M, Alshahni MM, Tanaka R, Yaguchi T, Bontems O, Salamin K, Fratti M, Monod M. 2017. Terbinafine Resistance of *Trichophyton* Clinical Isolates Caused by Specific Point Mutations in the Squalene Epoxidase Gene. *Antimicrob Agents Chemother* 61:10.1128/aac.00115-17. <https://doi.org/10.1128/aac.00115-17>
- Moreno-Sabater A, Normand A-C, Bidaud A-L, Cremer G, Foulet F, Brun S, Bonnal C, Aït-Ammar N, Jabet A, Ayachi A, et al. Terbinafine Resistance in Dermatophytes: A French Multicenter Prospective Study. *Journal of Fungi*. 2022; 8(3):220. <https://doi.org/10.3390/jof8030220>
- De Paepe, R., Normand, AC., Uhrlaß, S. *et al.* Resistance Profile, Terbinafine Resistance Screening and MALDI-TOF MS Identification of the Emerging Pathogen *Trichophyton indotineae*. *Mycopathologia* 189, 29 (2024). <https://doi.org/10.1007/s11046-024-00835-4>

Contact

Roelke De Paepe • roelke.depaepe@sciensano.be • +32 2 642 55 33